
#34118 HeiQ Pure – Worldwide*:





HeiQ required test to use this hangtag	QV: Initial and after 20x washes QC: Initial and after 10x washes ISO 20743 or ASTM E2149 or JIS L 1902 (Staphylococcus aureus): Initial: Log reduction >= 3 After washes: Log reduction >= 2.2 AATCC 100 (equation 2) Initial: >= 90% After washes: >= 85%
*For use in following countries	Worldwide, except for (including but not limited to): EU, EEA, UK, Switzerland and Canada
Languages	Multilingual (EN/FR/DE/ES/IT)
Used for following HeiQ products	HeiQ Pure SPQR HeiQ Pure TAG

A

It is the manufacturer's / brand's responsibility to fulfill the local labeling requirements and ensure that the correct hangtag is being attached to their product!